


```
In [21]: #Biomedical Signal Processing
```

```
In [7]: #import needed library
import numpy as np
import pandas as pd
import neurokit2 as nk
```

```
In [8]: # Generate synthetic signals
ecg = nk.ecg_simulate(duration=10, heart_rate=70)
ppg = nk.ppg_simulate(duration=10, heart_rate=70)
rsp = nk.rsp_simulate(duration=10, respiratory_rate=15)
eda = nk.eda_simulate(duration=10, scr_number=3)
emg = nk.emg_simulate(duration=10, burst_number=2)
```


```
In [9]: # Visualise biosignals
data = pd.DataFrame({"ECG": ecg,
                    "PPG": ppg,
                    "RSP": rsp,
                    "EDA": eda,
                    "EMG": emg})
nk.signal_plot(data, subplots=True)
```



```
In [12]: #Electrodermal Activity (EDA/GSR)
# Generate 10 seconds of EDA signal (recorded at 250 samples / second) with 2 SC
eda = nk.eda_simulate(duration=10, sampling_rate=250, scr_number=2, drift=0.01)

# Process it
signals, info = nk.eda_process(eda, sampling_rate=250)

# Visualise the processing
nk.eda_plot(signals, info)
```




```
In [13]: #Cardiac activity (ECG)
# Generate 15 seconds of ECG signal (recorded at 250 samples/second)
ecg = nk.ecg_simulate(duration=15, sampling_rate=250, heart_rate=70)

# Process it
signals, info = nk.ecg_process(ecg, sampling_rate=250)

# Visualise the processing
nk.ecg_plot(signals, info)
```

## Electrocardiogram (ECG)


ECG signal and individual Heart Beats (average heart rate: 69.9 bpm)



```
In [14]: #Respiration (RSP)
# Generate one minute of respiratory (RSP) signal (recorded at 250 samples / sec)
rsp = nk.rsp_simulate(duration=60, sampling_rate=250, respiratory_rate=15)

# Process it
signals, info = nk.rsp_process(rsp, sampling_rate=250)


# Visualise the processing
nk.rsp_plot(signals, info)
```



```
In [15]: #Photoplethysmography (PPG/BVP)
# Generate 15 seconds of PPG signal (recorded at 250 samples/second)
ppg = nk.ppg_simulate(duration=15, sampling_rate=250, heart_rate=70)

# Process it
signals, info = nk.ppg_process(ppg, sampling_rate=250)


# Visualize the processing
nk.ppg_plot(signals, info)
```



```
In [16]: #Electromyography (EMG)
# Generate 10 seconds of EMG signal (recorded at 250 samples/second)
emg = nk.emg_simulate(duration=10, sampling_rate=250, burst_number=3)

# Process it
signals, info = nk.emg_process(emg, sampling_rate=250)

# Visualise the processing
nk.emg_plot(signals, info)
```




```
In [20]: #Electrooculography (EOG)
# Import EOG data
eog_signal = nk.emg_simulate(duration=10, sampling_rate=250, burst_number=3)

# Process it
signals, info = nk.eog_process(eog_signal, sampling_rate=100)

# Plot
nk.eog_plot(signals, info)
```

# Electrooculography (EOG)

## Raw and Cleaned Signal



## Individual Blinks

